

2°C~8°C Refrigerator LT02-298U



LT02-298U



Specification:

Model	LT02-298U
External Size(W*D*H) mm	620*575*1980
Capacity	298L
Temp. Range	2°C~8°C
Temp. Accuracy	0.1°C
Control System	Microprocessor Control, Large Screen LED Display
Alarm	Audible and visual alarm for: High and low temperature, Power failure alarm, Door ajar, Thermostat failure
Refrigeration Type	Forced air refrigeration system
Refrigerant	R600a, CFC Free
Condenser & Evaporator	Bundy tube condenser
Defrost	Manual defrost / Auto defrost
Construction	Structure : Unibody Design & Mono-assembly foaming(Rigid Polyurethane insulation material) Internal : Embossed Aluminium Sheet External : Spraying Steel Plate
Shelves	5 pcs
Door	Glass door with electric heating
Consumption	189W
Power Supply	AC110/220V+10%, 50/60Hz
Standard Accessory	LED Lamp, USB port
Optional Accessory	Temperature recorder, Temperature printer(USB port can not be selected at the same time), probe access port
Package Size(W*D*H)mm	700*650*2030
Gross weight	84kg