

+2°C~+8°C

Medical Refrigerator



Constance

*Ultimate German Refrigeration
Technology*

LT02 Under-counter



Control System:

- Microprocessor controller, temperature ranging from +2°C to +8°C, can be set freely, controlling precision 1°C, display accuracy 0.1°C.
- Combined with two exact sensors and auto defrost.
- Audible and visual alarm: High or Low Temperature Alarm, Door Open Alarm, Sensor Failure Alarm.
- Standard with build-in USB datalogger, remote alarm contact and RS485 interface for monitor system.
- Power Supply: 220V / 50Hz, 1 phase.



Structure Design:

- Well designed in smart under-counter size for storage vaccine, biological products, reagent and drugs.
- Made of high quality steel, inner side with HIPS material, durable and easy to clean.
- With build-in LED lighting.
- Front opening lockable door with full height handle.
- With 2~3 High-quality steel wire shelves.
- Safety door lock, dual-glass-door design.
- 1 Cooling fan inside
- Optimal : Solid door



Refrigeration System:

- Germany EBM fan motor.
- Force air circulation system.
- With renowned brand compressor and condenser, better cool performance
- HCFC-FREE refrigerant ensure environmental protection & safety.



Specification:

Model	Capacity	External Size (WxDxH)mm	Input Power (W)	Weight (N.t) Kg
LT02-55U	55L	470 x 450 x 707	70	21
LT02-75U	75L	470 x 450 x 925	70	25