

-10°C~-25°C

Low Temperature Medical Freezer



Constance

*Ultimate German Refrigeration
Technology*



LT25 Series



Control System:

- Microprocessor temperature controller, temperature range: from -10°C to -25°C, can be set freely.
- Perfect audible/visual alarm: High or low temperature alarm, system failure alarm.
- Power Supply: 220V / 50Hz, 1 phase.



Structure Design:

- Upright type, ABS material interior, Outer body is painted steel board.
- 4pcs of caster.
- Safety lock on the door, with door handle.
- Drawer made of ABS materials



Refrigeration System:

- Germany EBM fan motor.
- R134a, CFC free refrigerant.
- Double compressor to control each chamber room separately (LT25-450U only)
- Certification: ISO9001, ISO14001, ISO13485.



Specification:

Model	Capacity	External Size (WxDxH)mm	Internal Size (WxDxH)mm	Input Power (W)	Weight (N.t) Kg
LT25-90U	90L	544 x 623 x 842	390 x 410 x 660	150	43
LT25-270U	270L	600 x 610 x 1750	410 x 414 x 1112	155	67
LT25-450U	450L	780 x 660 x 1850	640 x 515 x 1350	320	122