

-10°C~-25°C

Low Temperature Medical Freezer



Constance

*Ultimate German Refrigeration
Technology*

LT25 Chest Series



Control System:

- Microprocessor controller, temperature ranging from -10°C to -25°C, can be set freely, controlling precision 1°C, display accuracy 0.1°C.
- Delayed start and safe stop interval between restart and being terminated.
- Fast Freezing: Reach -25°C within 2-3 hours
- Audible and visual alarm: High or Low Temperature Alarm, Door Open Alarm, Sensor Failure Alarm.
- Power Supply: 220V / 50Hz, 1 phase.



Structure Design:

- Chest type, internal material are of stainless steel and outer are paint steel board.
- 4pcs of caster.
- Lockable door.
- Upward door opening with balancing door hinge
- 1 steel wire basket included.



Refrigeration System:

- R134a, CFC free refrigerant.
- Highly effective condenser and expansile evaporator to provide quick freezing.
- 40mm CFC free
- Certification: ISO9001, ISO14001, ISO13485.



Specification:

Model	Capacity	External Size (WxDxH)mm	Internal Size (WxDxH)mm	Input Power (W)	Weight (N.t/Gt) Kg
LT25-196C	196L	1035 x 630 x 928	870 x 400 x 510	71	39 / 50
LT25-358C	358L	1295 x 670 x 908	1080 x 525 x 728	121	56 / 70
LT25-508C	508L	1755 x 775 x 875	1570 x 620 x 640	210	75 / 90