

C350

2 in 1 Human Body and Industrial Thermometer



Main Feature

- High Temperature Alarm
- Temperature Difference can be Adjustable
- Laser Target Pointer Selection
- °C/°F Selection
- Backlight Display Selection
- Data Hold Function
- Auto Power Shut Off
- Low Battery Indication

Specification

1. Temperature Range
32°C~42°C
0°C~180°C
2. Accuracy: $\pm 0.3^{\circ}\text{C}/\pm 2\%$ or 2°C
3. Distance Spot Ratio: 5-15cm/12:1
4. Emissivity Adjustable: 0.95
5. Response Time & Wavelength: 0.8s & (8-14) μm
6. Repeatability: $\pm 1\%$ or $\pm 1^{\circ}\text{C}$
7. Resolution: 0.1°C or 0.1°F
8. Storage Temperature: $-20\text{--}50^{\circ}\text{C}$ ($-4\text{--}122^{\circ}\text{F}$)
9. Operating Temperature: $0\text{--}50^{\circ}\text{C}$ ($32\text{--}122^{\circ}\text{F}$)
10. Power supply: 1 pc DC 9V battery
11. Packaging Details : 20pcs/carton, carton size: 42*35*25cm, G.W.: 6kg